

Treating the World's Bumps & Bruises

RAPID
AID[®]

INFANT HEEL WARMER

A QUALITY HEEL STICK THE FIRST TIME



Our single use infant heel warmers are designed to provide the right amount of warmth to increase blood circulation in an infant heel during blood sampling. Made with our patented Insta-Gel technology, our heel warmers provide even heat distribution and consistent warming to help caregivers obtain a quality heel stick the first time.



BENEFITS

- **Effective:** Conforms to infant's foot shape, keeping warmth where needed
- **Safe:** Made with non-toxic food grade gel and boasts uniform heat distribution
- **Reliable:** Heel warmer will not exceed 40.5°C/105°F and gradually decreases temperature over time

FEATURES

- Comes in two easy activation methods:
 - Disk-activated
 - FTP™ (Fold-to-Pop) technology
- Delivers approximately 39.5°C/103°F of warm therapy when activated at room temperature
- Re-positionable adhesive strip holds warmer in place

✓ AVAILABLE AS PRIVATE LABEL OR
IN RAPID RELIEF BRAND

We have a variety of size and packaging options available to meet your needs.

For more information contact a Rapid Aid Account Manager:

Toll Free – 1-800-265-3468

Phone – 905-820-4788

Rapid Aid | 4120A Sladeview Crescent Units 1-4 Mississauga, ON L5L 5Z3

www.rapidaid.com

Find us on Facebook, LinkedIn, YouTube