

Treating the World's Bumps & Bruises



INSTANT COLD PACK

WITH GENTLE TOUCH™



Our Instant Cold Pack with Gentle Touch™ material provides on-the-spot cold therapy, featuring soft, non-woven material for direct-to-skin application. It works to immediately help reduce swelling and provide comfort after injury. More convenient than ice but just as effective, this portable pack can be used anywhere, anytime. Plus, with our FTP™ (Fold-to-Pop) technology, users can activate the pack simply by folding it in half.



BENEFITS

COLD THERAPY – For immediate injury treatment

- **Relieves aches & pains:** Reduces inflammation and muscle spasms
- **Reduces swelling & bruising:** Helps soothe pain by reducing blood flow to the area
- **Numbs the sore area:** Works effectively to reduce pain

FEATURES

- Gentle Touch™ material for direct-to-skin application
- Easy activation with FTP™ (Fold-to-Pop) technology
- Delivers approximately 2°C/36°F of cold therapy when activated at room temperature
- Choice of CAN, Urea, or Urea blend

✓ AVAILABLE AS PRIVATE LABEL OR IN RAPID RELIEF BRAND

We have a variety of size and packaging options available to meet your needs.

For more information contact a Rapid Aid Account Manager:
Toll Free – 1-800-265-3468
Phone – 905-820-4788

Rapid Aid | 4120A Sladeview Crescent Units 1-4 Mississauga, ON L5L 5Z3
www.rapidaid.com
Find us on Facebook, LinkedIn, YouTube