

RAPID RELIEF®
Stay in Motion

Medical Catalogue

Professional HOT & COLD
Therapy Products from
Concept to Delivery
– Branded or Private Label

CELEBRATING OVER 50 YEARS OF INNOVATION



Rapid Aid® is a privately owned company focused entirely on being the world leader and innovator in the temperature therapy market. We produce a wide array of branded and private label products for medical suppliers, retailers, and brands.

We cross multiple categories and utilize many technologies to deliver hot and cold therapy, to aid in recovery, relieve stress, and deliver pain relief.

Rapid Aid® can customize any one of our 500+ active skus or help you develop a completely new item that takes the industry by storm.

Rapid Aid® manages product development, testing, regulatory approval, package design, marketing assets, manufacturing, transportation, warehousing and distribution.

We have been in business over 50 years servicing the global market. With wholly owned manufacturing plants in Canada, US, China and Vietnam, you can rest easy knowing that inventory flow will be consistent no matter where your business is located.

Contact us today for consultation: info@rapidaid.com


The future of temperature therapy is allowing consumers and patients to safely, effectively and discreetly treat themselves hands free on-the-go. Just look for the Go-Therapy logo on available products.

GO THERAPY

Instant Cold Therapy	*4*
Instant Cold Packs	5
Instant Cold Packs with GENTLE TOUCH™	5
Instant Cold Packs Safety Kit	6
Disposable Ice Bags	7
Instant Hot Therapy	*8*
Instant Hot Packs	9
Instant Hot Packs with GENTLE TOUCH™	9
3-in-1 Instant Warm Wrap	9
Easy Draw Disc Activation ***NEW***	10
Easy Draw FTP Activation ***NEW***	10
Heat Bands	*11*
Air Activated Heat Bands	12
Air Activated Heat Bands	13
Cooling Gel Strips	*14*
Cooling Gel Strips Adult and Kids	15
Medicated Pain Relief Patches	*16*
4% Lidocaine Patch (U.S. only)	16
4% Lidocaine with 1% Menthol Patch (U.S. only)	16
5% Menthol Patch (U.S. only)	17
Plantar Fasciitis Medicated Patch (U.S. only)	17
5% Menthol Patch (Canada Only)	17
10% Methyl Salicylate Patch (Canada Only)	17
Reusable Therapy	*18*
Hot or Cold Compresses with Contour Gel	19
Hot or Cold Compresses with MAXFILL Contour Gel	19
Cold Vinyl Clay Gel Packs	20
Moist Heat Therapy Pads	21
Mom and Baby	*22*
Instant Cold Perineal Compress with Contoured Fit	23
Instant Cold Perineal Compress	23
Instant Warm Perineal Compress	23
Ultra-Thin Maxi Pad ***NEW***	24
Reusable Cool and Warm Gel Breast Compress	25
Reusable Cool Gel Breast Pads	25
Infant Heel Warmers	26
Infant Transport Mattress	26

RAPID AID Makes it Better!®

The broadest range of hot or cold therapy solutions in the industry, custom made to order



Wholly owned ISO facilities in Canada | USA | China | Vietnam

Multiple facility production ensures consistent inventory flow



INSTANT COLD THERAPY

Our Instant Cold Packs provide convenient, on-the-go cold therapy to help immediately reduce swelling and provide comfort after injury. Portable and easy to activate, they can be used anytime you need relief without the mess of ice cubes.

Available in Calcium Ammonium Nitrate (CAN) and Urea Blend (UB), they are designed to meet regional requirements while ensuring consistent performance.

Our Gentle Touch™ version features a soft outer layer, allowing for safe, direct-to-skin application without the need for a barrier, delivering the same effective relief wherever it is needed.

COLD THERAPY:
For immediate pain relief

- Relieves aches & pains: Reduces inflammation and muscle spasms
- Reduces swelling & bruising: Helps soothe pain by reducing blood flow to the area
- Numbs the sore area: Works effectively to reduce pain



Instant Cold Packs

Classic | Direct-to-Skin with Gentle Touch™







Item Number:	INC3124601PB /NA | INC0002201PB /EU	INC3125901PB / NA | INC0001601PB /NA	INC0002501PB /NA | INC0001701PB /EU
Size:	4.8″ x 6″	5.8″ x 8.5″	8″ x 11″
Material:	PE/PA	PE/PA	PE/PA
Fill:	Calcium Ammonium Nitrate	Calcium Ammonium Nitrate	Calcium Ammonium Nitrate

Classic Instant Cold Compress

- Delivers approximately 1˚C/34˚F of cold therapy when activated at room temperature
- Relieves aches & pains: Reduces inflammation and muscle spasms
- Reduces swelling & bruising: Helps soothe pain by reducing blood flow to the area
- Numbs the sore area: Works effectively to reduce pain







Item Number:	INC3185901PB /NA | INC0002401PB /EU	INC3135901PB /NA | INC0002101PB /NA	INC3134601PB /NA | INC0002001PB /EU
Size:	8.5″ x 5.8″	5.8″ x 8.5″	4.8″ x 6″
Material:	PE/PA/Non-Woven	PE/PA/Non-Woven	PE/PA/Non-Woven
Fill:	Calcium Ammonium Nitrate	Calcium Ammonium Nitrate	Calcium Ammonium Nitrate

Instant Cold Compress with Attached Strap

- All the same benefits of our Instant Cold Compress with Gentle Touch™, with the added benefit of a patented integrated strap.
- Conveniently acts as an instant cold pack, skin barrier and compression band all in one to provide all components of R.I.C.E. (Rest. Ice. Compression. Elevation.) therapy.

Instant Cold Compress with Gentle Touch™

- Delivers approximately 2˚C/36˚F of cold therapy when activated at room temperature.
- Provides all of the same benefits as our Classic version, but with our Gentle Touch™ material for safe, direct-to-skin application.



Instant Cold Pack

Safety Kits

Designed specifically for safety and first aid kits, our Instant Cold Pack is the most dependable choice for emergency cold therapy. It works to help immediately reduce swelling and provide comfort after injury. More convenient than ice but just as effective, this portable pack can be used anywhere, anytime.

GO THERAPY







Item Number:	INC0000501CB /NA | INC0001801CB /EU	INC0000601CB /NA | INC0001901CB /EU
Size:	4.8" x 6" (Product), 4" x 1.4" x 2.3" (Box)	5.8" x 8.5" (Product), 5" x 1.7" x 3.5" (Box)
Material:	PE/PA	PE/PA
Fill:	Calcium Ammonium Nitrate	Calcium Ammonium Nitrate

Temperature Performance – Rapid Relief Instant Cold Packs

This graph shows the average holding temperature of Rapid Relief Instant Cold Packs over a period of time, achieving an optimal temperature range for approximately 20 minutes.





Refillable Ice Bags

Fast relief for pain, swelling, bruising, and muscle cramps

These refillable bags feature a non-woven outer layer for comfort and reduced condensation. Free from natural rubber latex, they are safe for sensitive users. With sturdy ties and easy-to-use clamp closures, they are simple to secure and refill. Available in various sizes and designs to meet a range of patient needs.









Item Number:	IBC1000201PB	IBC1000301PB	IBC1000701PB	IBC1000501PB
Size:	7.8″ x 11.1″	5.2″ x 13.7″	4.6″ x 13.6″	5.2″ x 13.7″
Material:	PE/Non-Woven	PE/Non-Woven	PE/Non-Woven	PE/Non-Woven

Ice bag with Wire Closure For Eyes, Ears, Nose, Throat

- Ideal for patients recovering from eye or nasal surgery.
- Flexible wire closure for comfort.
- Universal size.

Ice bag with Clamp Closure

- Easy open-and-close reusable clamp closure.
- Universal size.


Perineal Ice Bag with Flexible Wire Closure

- Contoured shape for a good fit.
- Flexible wire closure for comfort.
- Universal size.

Dual Ice bag with Clamp Closure

- Bilateral design for treatment of maxillofacial and dental procedures.
- Easy open-and-close reusable clamp closure.
- Universal size.









Item Number:	IBC1000801PB	IBC1000401PB	IBC1000601PB	IBC1000901PB
Size:	7.3″ x 15.9″	5.2″ x 13.7″	5.2″ x 13.7″	7″ x 16.5″
Material:	PE/Non-Woven	PE/Non-Woven	PE/Non-Woven	PE/Non-Woven

Ice Bag with Clamp Closure, 4 Ties

- Four 14″ (35.6cm) ties allow bag to be attached right to affected area.
- Easy open-and-close reusable clamp closure
- Available in two sizes.

Ice Bag with Clamp Closure, Hook & Loop Straps

- Two hook & loop straps allow all bags to be attached right to the affected area.
- Easy open-and-close reusable clamp closure.
- Available in two sizes.

INSTANT HOT THERAPY


Our Instant Hot and Warm Packs, filled with patented Insta-Gel technology, deliver soothing heat instantly to promote therapeutic healing. Upon activation, the liquid turns into a contouring gel that molds to any area needing relief. While our Classic version is safe for direct-to-skin application, our Gentle Touch™ version features a soft, non-woven material that feels gentle on the skin, eliminating the need for a towel while providing comforting warmth where it's needed most.

HOT THERAPY RELIEF:

- Sore and stiff joints: Encourages relaxation and blood circulation
- Muscle aches and pain: Relaxes sore muscles by providing warm comfort
- Muscle tension and cramps: Increases blood flow to the affected area

Temperature Performance – Rapid Relief Instant Hot Pack

This graph shows the average holding temperature of Rapid Relief Instant Hot Packs over a period of time, achieving an optimal temperature range for approximately 23 minutes.





Instant Hot Packs

Classic | Direct-to-Skin with Gentle Touch™

Features

- Patented Insta-Gel technology turns from a liquid to contouring gel that molds to area being treated.
- Made with sodium acetate, it is designed to hit a specific temperature so heat can be controlled and customized.







Item Number:	INH4406601PB	INH4406901PB	INH4404901PB
Size:	5.8" x 5.8"	6" x 8.5"	4" x 9.3"
Material:	PE/PA	PE/PA	PE/PA
Fill:	Sodium Acetate & Water	Sodium Acetate & Water	Sodium Acetate & Water

Instant Hot Compress

- Delivers approximately 45˚C/113˚F of hot therapy when activated at room temperature.







Item Number:	INH4485901PB	INH4435901PB	INH4434601PB
Size:	5.8" x 8.5"	5.8" x 8.5"	4.8" x 6"
Material:	PE/PA/Non-Woven	PE/PA/Non-Woven	PE/PA/Non-Woven
Fill:	Sodium Acetate & Water	Sodium Acetate & Water	Sodium Acetate & Water

Instant Hot Compress with Attached Strap

- Warm pack with strap features built-in elastic band which allows for secure application and controlled compression (no need for bands or clips).

Instant Hot Compress with Gentle Touch™

- Delivers approximately 44˚C/111˚F of warm therapy when activated at room temperature.
- Gentle Touch™ material for safe, direct-to-skin application.



Easy Draw

Instant / Portable Hot Pack

Easy Draw™ Instant Hot Pack delivers fast, controlled heat on the go to help support blood draw and IV procedures. By promoting increased circulation and natural vasodilation at the application site, **Easy Draw™** helps make veins more visible and easier to access—reducing patient discomfort and improving procedural efficiency. Available in both **Fold-to-Pop** and **Disc-Activation** formats, **Easy Draw™** provides a reliable solution for clinical and home care environments.

Features:

- Delivers fast, controlled warmth using Fold-to-Pop or Disc-Activation technology for immediate use.

Benefits:

- Supports increased circulation at the application site
- Enhances vein visibility and palpability, helping reduce time spent searching for difficult veins.
- Promotes blood flow for capillary draws.
- Improves procedural efficiency while minimizing patient discomfort and multiple attempts.
- Enables smoother blood draws, reduces stress, and supports a more confident clinical workflow.
- Convenient, disposable, and travel-friendly.





Item Number	INH0001301PB (Fold to Activate)	INH0001201PB (Disc Activation)
Size	3.8″ x 5.5″	3.8″ x 5.5″
Material	PE/PA	PE/PA
Fill	Sodium Acetate & Water	Sodium Acetate & Water

EASIER VEIN IDENTIFICATION

Enhances vein visibility and palpability for both clinical and home use.

SUPPORTS FIRST-STICK SUCCESS

Improves procedural efficiency while minimizing patient discomfort and multiple attempts.

BOOST WORKFLOW & PATIENT EXPERIENCE

Enables easy blood draws. Ideal for in-home use kits.

HEAT BANDS

Ultra-thin, odorless Air Activated Heat Bands provide continuous, low-level heat for up to 8 hours, helping to relieve pain and soothe sore muscles. Designed for discreet wear under clothing, these non-medicinal, iron-filled wraps distribute heat evenly and are customizable in size, shape, color, and application to fit specific body areas.





Air Activated Heat Bands

Up to 8 Hours of hot therapy

Features

- Delivers long lasting heat therapy for up to 8 hours.
- Customized to fit the body area being treated.
- Ultra thin and can be worn discreetly under clothing.

Hot Therapy — For ongoing pain & stiffness

- Encourages quick recovery: Increases blood flow to help heal damaged tissue faster.
- Soothes stiff joints: Improves flexibility and helps improve mobility.
- Eases muscle aches: Adds comfort to muscles and calms soreness.

GO THERAPY





Item Number:	HBH1500401HP	HBH1500601HP
Size:	3.8" x 5.1"	3.8" x 10.3"
Material:	PE/Non-Woven	PE/Non-Woven
Fill:	Iron Filings	Iron Filings

Heat Band (Small)

- Compact size ideal for targeting small areas.
- Adheres directly to skin for hands-free application.

Heat Band (Large)

- Double the size of the small version, offering extended coverage for larger areas.
- Delivers long-lasting, even heat distribution for effective pain relief over a broader surface area.
- Adheres directly to skin for hands-free application.





Item Number:	HBH1500201HP	95303
Size:	L/XL	3.5" x 11.5"
Material:	PE/Non-Woven	PE/Non-Woven
Fill:	Iron Filings	Iron Filings

Heat Band Back & Hip

- Large, ergonomic design to cover both the lower back and hip regions.
- Ideal for providing continuous heat therapy to alleviate chronic back pain or hip discomfort.
- All-in-one wrap belt featuring hook & loop fastener to hold wrap in place.

Heat Band Neck & Shoulder

- Contoured design tailored to fit comfortably around the neck and shoulders.
- Provides targeted heat therapy to soothe muscle tension and stiffness in these common problem areas.
- Adheres directly to skin for hands-free application.



Air Activated Heat Bands

Up to 8 Hours of hot therapy



Item Number:	Neck & Shoulder Wrap
Size:	11.5″ x 3.5″
Material:	PE/Non-Woven
Fill:	Iron Filings

Neck & Shoulder Heat Band

- Contoured design tailored to fit comfortably around the neck and shoulders.
- Provides targeted heat therapy to soothe muscle tension and stiffness in these common problem areas.
- Adheres directly to skin for hands-free application.

Available for pre-order



Item Number:	Menstrual Heat Band
Size:	6.8″ x 3.7″
Material:	PE/Non-Woven
Fill:	Iron Filings

Heat Band Menstrual Pain Relief

- Air-activated warmth with iron fillings, no batteries or electricity required.
- Continuous heat therapy for up to 8 hours.
- Targeted relief for menstrual pain and muscle tension.

Available for pre-order



Item Number:	Menstrual Heart Shaped
Size:	3.9″ x 5.1″
Material:	PE/Non-Woven
Fill:	Iron Filings

Air Activated Menstrual Pain Relief Heat Band: Heart

- Air-activated warmth with iron fillings, no batteries or electricity required.
- Continuous heat therapy for up to 8 hours.
- Targeted relief for menstrual pain and muscle tension.
- Adheres to underwear for hands-free, secure application.
- Ultra-thin and discreet, ideal for wearing under clothing.
- Non-medicinal and fragrance-free, perfect for sensitive skin.

Available for pre-order

COOLING GEL STRIPS

Our cooling Gel Strips are medication free and safe to use with other prescriptions. Our unique combination of long-lasting cooling effect assists in relief from fever, migraines and hot flashes.

Pick your size and count to fit your needs.





Cooling Gel Strips

Direct-to-skin, non-medicated fever or migraine relief

Our Cooling Gel Strips provide cooling relief from fever, migraines, or hot flashes for up to 8 hours. Also effective in cooling the body by placing on the back of the neck, wrists, or under arms. These patches can be worn under clothing and can be easily applied directly to the skin adhering firmly, without leaving any sticky, greasy residue behind after removal. They are extremely flexible, soft and comfortable to wear.

GO THERAPY


Features

- Temperature: Delivers approximately 8 hours of cooling relief when used at room temperature.
- Ready to use.
- Direct-to-skin application.
- Peel-and-stick.
- No messy residue.

Cool Therapy Relief:


- Provides discreet cooling for up to 8 hours.
- Non-Medicated.
- Safe to use with medication.





Item Number:	96224 Kids 3ct	GSC2000103HP Adults 3ct
Size:	2" x 4.3"	2" x 5.1"
Material:	Gel + Non-Woven	Gel + Non-Woven
Fill:	Hydrogel	Hydrogel





MEDICATED PAIN RELIEF PATCHES

Our medicated patches offer targeted, penetrating therapy to relax muscles and soothe minor injuries. They are portable, convenient, and ready to use with both cooling and warmth effects. Easily applied directly to the skin with a peel-and-stick design, the patches adhere firmly without leaving sticky or greasy residue and can be worn discreetly under clothing.

Methyl Salicylate patches increase blood flow to promote faster recovery of damaged tissue.

Menthol patches provide localized relief by delivering menthol through the skin.

Lidocaine patches block pain signals at the nerve endings for effective local anesthesia.





Item Number:	MPN2500006CB	MPN2500406CB
Size:	4" x 5.5"	4" x 5.5"
Material:	Polypropylene / Non-Woven	Polypropylene / Non-Woven
Fill:	Hydrogel	Hydrogel

4% Lidocaine Medicated Patch

- Provides localized pain relief by numbing the nerves so they are less sensitive to pain.
- Active medical ingredient: 4% Lidocaine.

4% Lidocaine + 1% Menthol Medicated Patch


- Provides localized pain relief by numbing the nerves so they are less sensitive to pain.
- Added benefit of cooling sensation.
- Active medical ingredients: 4% Lidocaine, 1% Menthol.

Medicated Patches

Lidocaine, Menthol, Methyl Salicylate, Plantar Fasciitis







Item Number:	MPN2500808CB	MPC2500605CB
Size:	2.5" x 4"	3.9" x 7.9"
Material:	Polypropylene/Non Woven	Polypropylene/Non Woven
Fill:	Hydrogel	Hydrogel


Plantar Fasciitis

- Provides localized pain relief by numbing the nerves so they are less sensitive to pain.
- Active medical ingredient: 4% Lidocaine and 1% Menthol.

5% Menthol Medicated Patch

- Relieves aches & pains of the muscles/joints (such as arthritis, backache, and sprains).
- Numbs the sore area: Works effectively to reduce pain.
- Active medical ingredient 5% Menthol.





Item Number:	MPC2501306CB	MPH2502206CB
Size:	4" x 5.5"	4" x 5.5"
Material:	Polypropylene / Non-Woven	Polypropylene / Non-Woven
Fill:	Hydrogel	Hydrogel

5% Menthol Medicated Patch

- Relieves aches & pains of the muscles/joints (such as arthritis, backache, and sprains).
- Numbs the sore area: Works effectively to reduce pain.
- Active medical ingredient: 5% Menthol.

10% Methyl Salicylate Medicated Patch Encourages quick recovery: Increases blood flow to help heal damaged tissue faster.

- Soothes stiff joints: Improves flexibility and helps improve mobility.
- Eases muscle aches: Adds comfort to muscles and calms soreness.
- Active medical ingredient: 10% Methyl Salicylate.

REUSABLE THERAPY

Our reusable hot and cold therapy products are the most trusted and effective on the market. We provide an industry-leading selection of designs —from our classic reusable gel pack to our patented direct to skin premium fabric with inner liner to prevent leakage and extend product longevity.

COLD THERAPY

- Relieves aches & pains: Reduces inflammation and muscle spasms
- Reduces swelling & bruising: Helps soothe pain by reducing blood flow to the area
- Numbs the sore area: Works effectively to reduce pain

HOT THERAPY

- Encourages quick recovery: Increases blood flow to help heal damaged tissue faster
- Soothes stiff joints: Improves flexibility and help improve mobility
- Eases muscle aches: Adds comfort to muscles and calms soreness





Reusable Hot & Cold Compress

Contour-Gel | MAXFILL Contour-Gel

Our Reusable Hot & Cold Compresses are a must-have in every medicine cabinet, offering an economical solution for treating sports injuries, headaches, and general aches and pains. Our Contour-Gel molds to the body and remains flexible, when frozen, providing fast relief for inflammation. For enhanced comfort, the MAXFILL Contour-Gel Compress features added gel for extra weight, offering natural compression and extended relief. The soft, fully printed material ensures a quality feel, making it both effective and soothing while remaining flexibility even when frozen.

Features

- Delivers approximately -2°C/28°F of cold therapy when stored in the freezer (slightly lower or higher temperature depending on freezer time) and 43.5°C/110° F after microwaved.
- Different sizes available to use just about anywhere on the body from head to toe.
- Maintains flexibility when hot or cold.







Item Number:	RED1224601PB	RED1225901PB	RED1229001PB
Size:	4.8″ x 6″	6.3″ x 9″	9″ x 11″
Material:	PE/PA	PE/PA	PE/PA
Fill:	4 oz	10 oz	20 oz







Item Number:	RED0040101PB	RED0040201PB
Size:	6.4″ x 5″	10.6″ x 6″
Material:	PE/PA	PE/PA
Fill:	5 oz	12 oz





Item Number:	RED0040301PB	RED0040401PB	RED0040701PB
Size:	8.5″ x 6.8″	11.6″ x 10.8″	4.5″ x 23″
Material:	PE/PA	PE/PA	PE/PA
Fill:	12 oz	38 oz	19 oz



Cold Vinyl Clay Gel Packs

Our Cold Vinyl Clay Gel Packs are a staple in the medical and physiotherapy markets. An effective method for applying cold therapy to treat targeted areas of pain, vinyl clay gel packs remain flexible when frozen and wrap easily to the body part being treated. Easy to use, the clay clay gel is chilled in the freezer and maintains an optimal therapeutic temperature for up to 30 minutes.

Features

- Remains flexible when frozen
- Heavy duty and durable
- Made with non-toxic gel, natural clay
- Multiple shapes and sizes available to provide optimal coverage to the affected body areas being treated
- Versatile and flexible, conforms easily to the body for optimal cooling
- Clay gel evenly distributes cold therapy to soothe muscle tension and reduce inflammation







Item Number:	REC0502501PB	REC0502601PB	REC0502701PB
Size:	7.5" x 23"	6.5" x 7"	7" x 11.2"
Material:	PVC	PVC	PVC
Fill:	Clay-Based Gel	Clay-Based Gel	Clay-Based Gel
	Neck & Shoulders Cold Vinyl Clay Gel Wrap	**Quarter Size Cold Vinyl Clay Gel Wrap**	**Half Size Cold Vinyl Clay Gel Wrap**







Item Number:	REC0502201PB	REC0502301PB	REC0502401PB
Size:	11" x 14"	11" x 21"	4" x 12"
Material:	PVC	PVC	PVC
Fill:	Clay-Based Gel	Clay-Based Gel	Clay-Based Gel
	Standard Cold Vinyl Clay Gel Wrap	**Oversize Cold Vinyl Clay Gel Wrap**	**Slim Cold Vinyl Clay Gel Wrap**



Moist Heat Therapy Pads

Trusted by medical professionals for decades, Moist Heat Pads provide penetrating moist heat to relieve tension and soothe sore muscles, and provide pain relief due to arthritis, back pain, muscle strains and more.

Features

- Provides up to 30 minutes of deep, soothing moist heat
- Different colored tabs for easy rotation in busy clinics and easy removal from the heating unit
- All-natural filled pack wrapped in a woven cotton blend
- Variety of sizes to accommodate different body parts
- Offers deep, penetrating heat to alleviate pain and promote relaxation
- Soothing heat therapy to relieve tension and stiffness





Item Number:	MHH4500001PB	MHH4500101PB
Size:	5.1″ x 11.8″	7.1″ x 23.6″
Material:	Cotton Blend	Cotton Blend
Fill:	Silica Clay Granules	Silica Clay Granules
	Half Size Hydrocollator Moist Heat Pack	**Neck & Shoulders Hydrocollator Moist Heat Pack**







Item Number:	MHH4500201PB	MHH4500301PB	MHH4500401PB
Size:	9.8″ x 20.1″	9.8″ x 11.8″	15″ x 24″
Material:	Cotton Blend	Cotton Blend	Cotton Blend
Fill:	Silica Clay Granules	Silica Clay Granules	Silica Clay Granules
	Knee & Shoulder Hydrocollator Moist Heat Pack	**Standard Size Hydrocollator Moist Heat Pack**	**Oversize Hydrocollator Moist Heat Pack**



Mom & Baby

The Instant Cold Perineal Compress (direct-to-skin) delivers quick cold therapy to ease pain and reduce swelling in the perineal area, offering maximum comfort with minimal nursing intervention. The hourglass shape of our Premium Perineal Compress contours to the body while providing enhanced absorbency.

The Instant Warm Perineal Compress helps soothe perineal pain before and after childbirth by providing comforting warmth, promoting tissue relaxation, and improving circulation. It is a convenient alternative to a sitz bath with similar effectiveness.

Our non-toxic, phthalate-free, BPA-free, and reusable gel pads and compresses are soft, flexible, and contoured to comfortably mold to the breast. With extra soft covers, they can be worn discreetly under a bra.

Our patented Infant Heel Warmers provides consistent temperature immediately upon activation and contours to the baby's heel to increase blood flow, resulting in a more comfortable experience for the infant and more accurate heel stick conducted by caregivers.

Made with our patented Insta-Gel, our mattress helps bring warm comfort to infants during transport in or between hospitals, reducing the risk of cold stress to infants.



Instant Perineal Compress

Direct to Skin

Features

- Extra absorbent pad with a long front-to-back area for added coverage and an adhesive strip to keep it firmly in place
- FTP (Fold-to-Pop) technology for easy activation
- Not made with natural rubber latex
- Direct-to-skin application

Instant Cold Perineal Compress provides fast-acting cold therapy to relieve pain and reduce swelling in the perineal area with minimal nursing intervention, making it ideal for postpartum care and post-vasectomy recovery.







- **Extra wide contour shape**
- **Optimal absorbancy**
- **No need for extra OB pads**

Item Number:	PPC9164001PB /NA | PPC3500401PB /EU	PPC9166001PB /NA | PPC3500501PB /EU
Size:	4.6" x 13.8"	7.5" x 14.5"
Material:	Fluff Pulp/Non-Woven	Fluff Pulp/Non-Woven
Fill:	Calcium Ammonium Nitrate	Calcium Ammonium Nitrate

Instant Cold Perineal Compress

- Delivers approximately 4.4˚C/40˚F of cold therapy when activated at room temperature.
- Relief from pain and inflammation of the perineal area after childbirth.
- Effective also for post vasectomy procedures.
- Unique hourglass shape of premium perineal pad is designed to contour the body's natural curves while delivering more absorbency.

Instant Warm Perineal Compress delivers soothing warmth to reduce perineal pain before and after childbirth, promoting tissue relaxation and improved circulation. Convenient and effective, it is an excellent alternative to a warm sitz bath and can also help relieve menstrual cramps.



Item Number:	PPH9464001PB
Size:	4.6" x 13.8"
Material:	Fluff Pulp/Non-Woven
Fill:	Sodium Acetate & Water

Instant Warm Perineal Compress

- Delivers approximately 45.6˚C/114˚F of warm therapy when activated at room temperature.
- Soothes perineal area.
- Encourages relaxation and blood circulation prior to childbirth (typically 72 hours warm therapy use)
- Relieves tension: Can be used to soothe menstrual cramps.

Instant Cold Ultra Thin Maxi-Pad

Direct to Skin

Patent Pending 19/106,516




Our Ultra Thin Instant Cold Maxi Pad provides the same fast-acting relief and leakage protection from perineal pain after childbirth. This Ultra Thin Maxi-Pad is designed for comfort first and can be use direct-to-skin. It features an extra-absorbent pad that is 26% thinner and 30% lighter than the leading national brand. Easy activation via FTP™ technology, and controlled cold therapy at approximately -7˚C/20˚F to 8˚C/47˚F. Latex-free and available in bi-lingual flexible outer packaging.

Features:

- Ultra Thin. 26% Thinner, and 30% lighter than the leading national brand.
- Delivers between -7˚C/20˚F to 8˚C/47˚F of cold therapy when activated at room temperature.
- Direct-to-skin application.
- Extra absorbent pad with an adhesive strip to keep it firmly in place.
- FTP™ (Fold-to-Pop) technology for easy activation.
- Not made with natural rubber latex.

Fast relief for:

- Pain and inflammation of the perineal area after childbirth.
- Effective also for post vasectomy procedures.







Comfort First - Ultra Thin

- 26% thinner and 30% lighter than the leading brand.

Item Number:	PPC3500308CB EU | PPC3501008CB NA
Size:	4.6" x 11.5"
Material:	Fluff Pulp/Non-Woven
Fill:	Calcium Ammonium Nitrate

Instant Cold Ultra-Thin Perineal Maxi-pad

- Reliable absorption without the additional bulk.
- Soothing cooling center strip for increased contact where needed most.
- Functional, slimmer design allows for enhanced comfort.
- Compact design delivers smaller shelf foot-print.
- Reduced size provides significant logistical cost savings.



Reusable Cool & Warm Gel Breast Compress

Our Reusable Cool Gel Breast Pads reduce discomfort before and after a biopsy, mastectomy, or cosmetic surgery. Their slim design allows discreet wear, while the cool gel naturally conforms to the patient's anatomy.

Our Cool and Warm Gel Breast Compresses provide fast relief through cool or warm therapy, helping to reduce swelling, ease pain during breastfeeding or post-surgery, and support lactation.

Features

- Contour-Gel technology allows the pad to gently contour to the body area being treated and remain flexible when frozen
- Direct-to-skin application with Gentle Touch™ technology







Item Number:	FCC3000701PB	FCC3000601PB	FCD3000501PB
Size:	3" x 3"	5.3" x 4.8"	7" x 7"
Material:	PE/PA/Non-Woven	PE/PA/Non-Woven	PE/PA/Non-Woven
Fill:	Contour-Gel	Contour-Gel	Contour-Gel

Reusable Cool Gel Breast Pads

- Delivers approximately 1°C/34°F of cool relief when stored in the freezer (slightly lower or higher temperature depending on freezer time)
- Slim design can fit inside the bra

Cool Therapy Relief

- Relieves pain due to post-cosmetic surgical swelling
- Helps alleviate pain and discomfort post-mastectomy
- Minimizes post-biopsy swelling and bruising

Reusable Cool & Warm Gel Breast Compress

- Delivers approximately 42°C/108°F of warm therapy when heated in a microwave, and approximately 1°C/34°F of cool relief when stored in the freezer (slightly lower or higher temperature depending on freezer time). Temperature adjustable based on user's preference.

Cool Therapy Relief

- Relieves pain and swelling due to engorgement
- Helps alleviate blocked ducts and mastitis pain

Warm Therapy Relief

- Relieves engorgement, plugged ducts, mastitis or to encourage let-down
- Use with breast pump to encourage let-down and reduce the time spent pumping



Infant Care

Heel Warmers | Transport Mattress

Our Infant Heel Warmer features patented Insta-Gel technology, activating instantly to provide consistent warmth and contour to the baby's heel. With direct-to-skin application, it improves blood flow, making vein access easier and creating a more comfortable experience for infants during heel sticks.

Our gel-based Infant Transport Mattress ensures warm comfort for infants during hospital or EMS transport. Upon activation, it maintains a stable temperature, reducing the risk of cold stress, while its Insta-Gel design provides safety and comfort when needed most.







Item Number:	IFW9423501PB (Fold Activation)	IFW9424501PB (Disc Activation)	IFW9165001PB (Disc Activation)
Size:	3.8" x 5.5"	3.8" x 5.5"	9.7" x 15.8"
Material:	PE/PA	PE/PA	PE/PA/Non-Woven
Fill:	Sodium Acetate & Water	Sodium Acetate & Water	Sodium Acetate and Water

Infant Heel Warmer

- Delivers approximately 39.5°C/103°F heat when activated at room temperature
- Adhesive strip that holds the product in place
- Fold-to-Pop technology, allowing the user to activate the heel warmer simply by folding it in half
- Not made with natural rubber latex
- Patented formula prevents pre-activation during transport

Infant Heel Warmer

- Delivers approximately 39.5°C/103°F heat when activated at room temperature
- Adhesive strip that holds the product in place
- Easy activation disc
- Not made with natural rubber latex

Infant Transport Mattress

- Delivers approximately 39.5°C/103°F of warm therapy when activated at room temperature and will not exceed 40.5°C/104°F
- Single-use minimizes the risk of infection
- Portable heat source, no electricity required
- Easy activation disc
- Not made with natural rubber latex

RAPID RELIEF®
Stay in Motion

4120A Sladeview Crescent
Mississauga, ON Canada L5L 5Z3
www.rapidaid.com